Medical image enhancement

For whom, why, and how to.

For whom?

- Human processing
- Computer processing

Why?

- Can't distinguish between tissues
- Data too noisy for computer algorithm to perform well
- Acquisition or reconstruction artifacts interfere with visualization or algorithm processing

2D/3D image processing tools

Many general purpose tools (or filters) exist to remove noise, increase contrast, etc.

Most IP applications must employ a suite of such tools to achieve their goal.

Linear, shift-invariant filters

- Obey principle of superposition F(x) + F(y) = F(x+y)
- Are not position-dependent
- Can be implemented in either spatial domain (by convolution) or frequency domain (by weighting components)

How to filter to

- Remove noise
- Emphasize edges
- Detect edges
- Modify shapes

Increase contrast

- Scaling (shift-invariant, linear)
- Window/level
- Histogram equalization (shift-invariant, but nonlinear)
- Adaptive histogram equalization

Increase contrast - scaling

Scaling simply spreads or contracts image intensity values into a new range – e.g. to monitor display range.


Scaling, if the final output is required to have integer values, may not be a one-to-one operation.

Before scaling, consider cropping the image to the region of interest. This can also help your visual system.

Increase contrast – linear rescaling

Digitized film, 12 bit storage, scaled to 256 colors






Increase contrast – linear rescaling

Simply-cropped image, 12 bit storage, scaled to 256 colors





Increase contrast — window/level

Window and level adjustments are piecewise-linear or nonlinear operations that can be approximated by adjusting the lookup table of a displayed image.


Except on specialized workstations, this will not recover contrast that is lost when the data range of an image was reduced for display.

More properly, contrast adjustment is performed on original data using a lookup table and then mapped to a display.

Increase contrast - window/level

Window is narrowed to span less of input data; level (midpoint of slope) is shifted right to decrease brightness.






Increase contrast - window/level


Window is narrowed to span less of input data; level (midpoint of slope) is shifted left to increase brightness.






Increase contrast – histograms

Histogram equalization attempts to take advantage of unassigned output values. It is a one-to-one operation, unlike window/level operations. The lookup table simply follows the shape of the cumulative histogram.






Increase contrast – histogram equalization

Raw image

Global HEQ





Increase contrast – adaptive histogram equalization

Adaptive histogram equalization computes an output intensity level from the histogram of a local neighborhood of each pixel. Many variants exist that attempt to control contrast, interpolate for speed, etc.

Global HEQ

Fully adaptive HEQ



Increase contrast – adaptive histogram equalization

Contrast limitation in histogram equalization means setting a maximum spread of output intensities for adjacent input intensities.

Raw image

Contrast-limited fully adaptive HEQ



How to filter to

- ∠ Increase contrast
- Remove noise
- Emphasize edges
- Detect edges
- Modify shapes


Simple convolution filters (linear)

- Replace center pixel with weighted combination of local neighbors; a moving average
- Mean or "boxcar" filter; uniform weights
- Gaussian kernel; center weights dominate


Test object

With additive noise


Boxcar mean, 5x5




Plot shows histograms of test object and filtered image. Additive noise is brighter than all objects.



Test object With additive noise Gaussian mean, fwhm 5



Distribution of intensity values after smoothing differs from simple mean.



Trimmed mean filter (nonlinear)


- Boxcar kernel
- Discard highest and lowest values
- Compute mean of remaining values

Sometimes called an "Olympic" filter


Test object

With additive noise

Trimmed mean, 5x5




Removes most noise, with less effect on primary signal intensity.





T1-weighted MRI slice






Distinct histogram peaks in original are poorly differentiated after smoothing.







T1-weighted MRI slice






Note change in histogram peaks after smoothing.







T1-weighted MRI slice









Median filters (non-linear)

- Gather N values using a kernel
- Sort by value
- Replace center pixel by median value

Removes features within the kernel that occupy less than N/2 pixels, also trims off edges of larger features. Can leave artifacts that reflect the kernel shape.

Median kernels; omni-directional

1	1	1	1	1
1	1	1	1	1
1	1	1	1	1
1	1	1	1	1
1	1	1	1	1

	1	1	1	
1	1	1	1	1
1	1	1	1	1
1	1	1	1	1
	1	1	1	


		1		
	1	1	1	
1	1	1	1	1
	1	1	1	
		1		

Center-weighting can be achieved by adding 2K copies of the central pixel prior to taking the median value.


Test object

With additive noise

Square median, 5x5




Noise is neatly suppressed – but also the open circle.




Test object

Center-weighted median, 5x5




Center pixel got 3 votes instead of 1. No discernable difference from simple median.




T1-weighted MRI slice






Note that bright rim, midhemisphere boundary have been suppressed.


Histogram peaks are accentuated, rather than smeared together.







T1-weighted MRI slice

Center-weighted median, 5x5



Only subtly differs from ordinary median.






Median kernels; directional


	1	
	1	
	1	
	1	
	1	

1		
	1	
	1	
	1	
1		




Test object

5-point vertical line median




Retains a bit of the open circle; small signals keep their shape; a few noise specks remain.




Test object

5-point curve median




Retains more of the open circle; a few more noise specks remain.




Test object

5-point rotated curve median




Retains different portions of the open circle; different noise specks survive.




Noise removal – median filters


T1-weighted MRI slice






Rim is less degraded, histogram is similar to original.






Noise removal – median filters


T1-weighted MRI slice



Similar to line kernel; left and right rims are well preserved.


5-point curve median





Noise removal – median filters


T1-weighted MRI slice




Now it's the fore and aft rim that is well-preserved.

Note that mid-hemisphere boundary is less distinct.

5-point rotated curve median





Noise removal — Wiener filter

Wiener filter attempts to minimize mean-square-error

$$g = f + n \qquad \hat{f} = h * g \qquad \hat{F} = HG$$

$$H = S_{ff}$$

$$S_{ff} + S_{nn}$$

where S_{ff} is the Fourier transform of the autocorrelation function of the image and, if we have stationary white noise, $S_{nn} = \frac{s^2}{n}$ (known or estimated from a section of the image)

Noise removal - anisotropic smoothing

2D iterative method for "diffusion-weighted" smoothing

	N	
Е	Т	W
	S	

$$T' = T + p*(Cn*(N-T) + Cs*(S-T)) + p*(Ce*(E-T) + Cw*(W-T))$$


$$Cn = exp(-1 * (|N-T|/K)^2)$$

When gradient is large, conduction coefficient is negligible, so that neighbor can not affect the value of T. Extend to 3D by including Up and Down terms.

Noise removal - anisotropic smoothing

Raw

Anisotropic smoothing, K=10, iterations = 6



How to filter to

- ∠ Increase contrast
- Remove noise
- Emphasize edges
- Detect edges
- Modify shapes


Emphasize edges — spatial domain

Edge-boost methods use directional or omnidirectional kernels to identify edges; these can be recombined as

g = image + a(edges)

Unsharp masking uses a convolution kernel to compute a very blurred version of an image; the enhanced image is defined as g=unsharp + b(image-unsharp)

Emphasize edges – spatial domain



Emphasize edges – frequency domain

Compute Fourier transform of image


Apply a weighting function to boost high frequency components

Compute inverse transform

Emphasize edges – frequency domain

Original, boosted

FFT (original, modified)



How to filter to

- ∠ Increase contrast
- Remove noise
- Emphasize edges
- Detect edges
- Modify shapes

- Horizontal, vertical, omni-directional kernels
- Sobel filter
- Difference of Gaussians
- Marr-Hildreth filter

Convolution kernels that take derivatives

Gx = horizontal derivative: Gy = vertical derivative:


-1	0	1
-2	0	2
-1	0	1

-1	-2	-1
0	0	0
1	2	1


Laplacian = second derivative

-1	-1	-1
-1	8	-1
-1	-1	-1


Edge detection – horizontal derivative



Edge detection – vertical derivative




Edge detection – Laplacian



Nonlinear methods often combine directional derivatives.

Roberts' cross:
$$_2$$
 result = SQRT([F - A4] + [A3-A5])

Edge detection – Roberts'





Pixel labels (2D):

A0	A1	A2
A7	F	A3
A6	A5	A4

Sobel operator: result = SQRT(X*X + Y*Y) where X = (A2+2A3+A4) - (A0+2A7+A6)Y = (A0+2A1+A2) - (A6+2A5+A4)

Edge detection – Sobel





Difference of Gaussians ("DOG")


result = Gauss(s1)**image - Gauss(s2)**image

This is believed to be similar to processing in the human

visual system.

Edge detection - DOG

Gaussian 1 Gaussian 2 Difference



Marr-Hildreth operator -- filter is a Laplacian of a Gaussian

2D:
$$C(x,y) = \nabla (I(x,y) *** G(x,y,r))$$

3D: $C(x,y,z) = \nabla (I(x,y,z) *** G(x,y,z,s))$
 $= I(x,y,z) *** \nabla (G(x,y,z,s))$

Result is very similar to DOG filtering

Edge detection – zero crossings

2D pixel labeling:

A0	A1	A2
A7	F	A3
A6	A5	A4


```
If A0 & F are different signs, Grad1 = |A0 - F|, else Grad1 = 0.0
If A1 & F are different signs, Grad2 = |A1 - F|, else Grad2 = 0.0
```

٠

Result = Maximum of Grad1 ... Grad8

Original

Filtered with 2D Marr-Hildreth operator
Zero-crossings
Zero-crossings overlaid on original data




Frei and Chen filter:

This sums the result of 9 convolution filters, then takes the COS(SQRT(sum)).

The first kernel is a boxcar average; the others are various derivative kernels

Ref. J.C. Russ, *The Image Processing Handbook*, CRC Press, 1995

Edge detection – Frei and Chen





How to filter to ...

- ∠ Increase contrast
- Remove noise
- Emphasize edges
- Detect edges
- Modify shapes

Shape modification tools

Morphological (shape-based) operations

- Erosion
- Dilation
- Opening
- Closing

Have predictable effects on shape.


Morphological operators

Erosion ("thin" or "remove outer layer") ⊗

- -- Gather N values using a kernel
- -- Sort by value
- -- Replace center pixel by minimum value

Islands smaller than the kernel are removed. Small pennisulas are broken.

Morphological operators - erosion



raw,
3x3 erosion


eroded 2 times, eroded 3 times

Morphological operators

- -- Gather N values using a kernel
- -- Sort by value
- -- Replace center pixel by maximum value

Holes smaller than the kernel are filled. Neighboring islands may be connected.

Morphological operators - dilation



raw, 3x3 dilation


dilated 2 times, dilated 3 times

Morphological operators

Opening = dilation followed by erosion

Smoothes a contour
Breaks narrow isthmuses
Eliminates small islands
Eliminates sharp capes

Morphological operators - opening



Morphological operators

Closing = erosion followed by dilation

Smoothes a contour Fuses narrow breaks Fills small holes Fills small gaps on a contour

Morphological operators - closing

